

IFT81 Monoclonal antibody, PBS Only

Catalog Number: 67987-1-PBS

Basic Information

Catalog Number:	67987-1-PBS	GenBank Accession Number:	BC029349	Purification Method:	Protein G purification
Size:	1 mg/ml	GeneID (NCBI):	28981	CloneNo.:	1D5A3
Source:	Mouse	UNIPROT ID:	Q8WYAO		
Isotype:	IgG1	Full Name:	intraflagellar transport 81 homolog (Chlamydomonas)		
Immunogen Catalog Number:	AG31629	Calculated MW:	676 aa, 80 kDa		
		Observed MW:	70-80 kDa		

Applications

Tested Applications:
WB, Indirect ELISA, IF

Species Specificity:
Human, Rat, Mouse, Rabbit, canine

Background Information

Intraflagellar transport (IFT), mediated by molecular motors and IFT particles, is an important transport process that occurs in the cilium and has been shown to be essential for the assembly and maintenance of cilia and flagella in many organisms. IFT particles are multi-subunit complexes of proteins that function to move non-membrane-bound particles from the cell body to the tip of cilium or flagellum, then return them to the cell body. Transport towards the ciliary tip is regulated by the IFT complex B (IFT-B), consisting of at least 15 IFT proteins, in association with kinesin motors, whereas transport from the ciliary tip back to the base is executed by a dynein motor in association with the IFT complex A (IFT-A), currently known to be composed of six IFT proteins. IFT81 is a subunit of IFT complex B. It may play a role in development of the testis and spermatogenesis. There are some isoforms of IFT81 with 73-78 kDa and 43-50 kDa.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 67987-1-Ig (IFT81 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67987-1-PBS in a different storage buffer formulation.



rabbit brain tissue were subjected to SDS PAGE followed by western blot with 67987-1-Ig (IFT81 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67987-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 67987-1-Ig (IFT81 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67987-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed hTERT-RPE1 cells using IFT81 antibody (67987-1-Ig, Clone: 1D5A3) at dilution of 1:2000 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L), ARL13B antibody (17711-1-AP, red). This data was developed using the same antibody clone with 67987-1-PBS in a different storage buffer formulation.