

For Research Use Only

TGM2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68006-2-PBS



Basic Information

Catalog Number:

68006-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7462

GenBank Accession Number:

BC003551

GeneID (NCBI):

7052

ENSEMBL Gene ID:

ENSG00000198959

UNIPROT ID:

P21980

Full Name:

transglutaminase 2 (C polypeptide, protein-glutamine-gamma-glutamyltransferase)

Calculated MW:

77 kDa

Purification Method:

Protein A purification

CloneNo.:

2B3G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51456-1, TGM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68006-1-PBS. Detection antibody: 68006-2-PBS.



Cytometric bead array standard curve of MP51456-1, TGM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68006-1-PBS. Detection antibody: 68006-2-PBS. Standard: Ag7462. Range: 3.125-400 ng/mL.