

For Research Use Only

# DNM1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68009-3-PBS



## Basic Information

Catalog Number:

68009-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12932

GenBank Accession Number:

BC063850

GeneID (NCBI):

1759

UNIPROT ID:

Q05193

Full Name:

dynamin 1

Calculated MW:

851 aa, 96 kDa

Purification Method:

Protein G purification

CloneNo.:

2D8C3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sample test of MP51318-1, DNM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68009-2-PBS. Detection antibody: 68009-3-PBS.



Cytometric bead array standard curve of MP51318-1, DNM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68009-2-PBS. Detection antibody: 68009-3-PBS. Standard:Ag12932. Range: 0.977-500 ng/mL.