

For Research Use Only

# DNM1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68009-4-PBS



## Basic Information

|                                      |                                       |  |
|--------------------------------------|---------------------------------------|--|
| Catalog Number:<br>68009-4-PBS       | GenBank Accession Number:<br>BC063850 | Purification Method:<br>Protein G purification |
| Concentration:<br>1 mg/ml            | GeneID (NCBI):<br>1759                | CloneNo.:<br>2A9H9                             |
| Source:<br>Mouse                     | UNIPROT ID:<br>Q05193                 |  |
| Isotype:<br>IgG1                     | Full Name:<br>dynamin 1               |  |
| Immunogen Catalog Number:<br>AG12932 | Calculated MW:<br>851 aa, 96 kDa      |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Sample test of MP51318-2, DNM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68009-2-PBS. Detection antibody: 68009-4-PBS.



Cytometric bead array standard curve of MP51318-2, DNM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68009-2-PBS. Detection antibody: 68009-4-PBS. Standard:Ag12932. Range: 0.977-500 ng/mL.