

For Research Use Only

RALY Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68011-2-PBS



Basic Information

Catalog Number: 68011-2-PBS	GenBank Accession Number: BC103753	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 22913	CloneNo.: 3A10D7
Source: Mouse	UNIPROT ID: Q9UKM9	
Isotype: IgG1	Full Name: RNA binding protein, autoantigenic (hnRNP-associated with lethal yellow homolog (mouse))	
Immunogen Catalog Number: AG12250	Calculated MW: 307 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

