

For Research Use Only

# RALY Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68011-3-PBS



## Basic Information

<b>Catalog Number:</b> 68011-3-PBS	<b>GenBank Accession Number:</b> BC103753	<b>Purification Method:</b> Protein A Magarose purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 22913	<b>CloneNo.:</b> 2H1H10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9UKM9	
<b>Isotype:</b> IgG3	<b>Full Name:</b> RNA binding protein, autoantigenic (hnRNP-associated with lethal yellow homolog (mouse))	
<b>Immunogen Catalog Number:</b> AG12250	<b>Calculated MW:</b> 307 aa, 33 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

