

For Research Use Only

GNAO1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68027-2-PBS



Basic Information

Catalog Number:

68027-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG32045

GenBank Accession Number:

BC030027

GeneID (NCBI):

2775

UNIPROT ID:

P09471

Full Name:

guanine nucleotide binding protein (G protein), alpha activating activity polypeptide O

Calculated MW:

302 aa, 35 kDa

Purification Method:

Protein G purification

CloneNo.:

1A1B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51895-1, GNAO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68027-1-PBS. Detection antibody: 68027-2-PBS.



Cytometric bead array standard curve of MP51895-1, GNAO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68027-1-PBS. Detection antibody: 68027-2-PBS. Standard: Ag32045. Range: 1.563-100 ng/mL.