

For Research Use Only

# Endophilin 1 Monoclonal antibody, PBS Only



Catalog Number: 68029-1-PBS

## Basic Information

Catalog Number:

68029-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG24793

GenBank Accession Number:

BC032825

GeneID (NCBI):

6456

UNIPROT ID:

Q99962

Full Name:

SH3-domain GRB2-like 2

Calculated MW:

39.8 kDa

Observed MW:

38-40 kDa

Purification Method:

Protein A purification

CloneNo.:

1G4H1

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat, pig, rabbit

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68029-1-Ig (SH3GL2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68029-1-PBS in a different storage buffer formulation.