

For Research Use Only

PAIP1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68042-2-PBS



Basic Information

Catalog Number: 68042-2-PBS	GenBank Accession Number: BC005295	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 10605	CloneNo.: 2D7B8
Isotype: IgG1	UNIPROT ID: Q9H074	
Immunogen Catalog Number: AG16773	Full Name: poly(A) binding protein interacting protein 1	
	Calculated MW: 70 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51721-1, PAIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68042-2-PBS. Detection antibody: 68042-3-PBS.



Cytometric bead array standard curve of MP51721-1, PAIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68042-2-PBS. Detection antibody: 68042-3-PBS. Standard: Ag16773. Range: 3.13-50 ng/mL