

For Research Use Only

TTLL12 Monoclonal antibody, PBS Only



Catalog Number: 68044-1-PBS

Basic Information

Catalog Number:

68044-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13804

GenBank Accession Number:

BC001070

GeneID (NCBI):

23170

UNIPROT ID:

Q14166

Full Name:

tubulin tyrosine ligase-like family, member 12

Calculated MW:

644 aa, 74 kDa

Observed MW:

72-74 kDa

Purification Method:

Protein G purification

CloneNo.:

1H5C11

Applications

Tested Applications:

WB, Indirect ELISA, IHC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68044-1-Ig (TTLL12 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68044-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 68044-1-Ig (TTLL12 antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 68044-1-PBS in a different storage buffer formulation.