

For Research Use Only

SH3BP4 Monoclonal antibody

Catalog Number: 68045-1-Ig



Basic Information

Catalog Number:

68045-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31238

GenBank Accession Number:

BC057396

GeneID (NCBI):

23677

UNIPROT ID:

Q9POV3

Full Name:

SH3-domain binding protein 4

Calculated MW:

963 aa, 108 kDa

Observed MW:

110 kDa

Purification Method:

Protein G purification

CloneNo.:

2G4H3

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : Y79 cells, U2OS cells, A549 cells, HeLa cells, HepG2 cells, hTERT-RPE1 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various cell lysates were subjected to SDS PAGE followed by western blot with 68045-1-ig (SH3BP4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with Lamin B1 Monoclonal antibody (66095-1-ig) as loading control.