

For Research Use Only

SH3BP4 Monoclonal antibody, PBS Only



Catalog Number: 68045-1-PBS

Basic Information

Catalog Number: 68045-1-PBS	GenBank Accession Number: BC057396	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 23677	CloneNo.: 2G4H3
Source: Mouse	UNIPROT ID: Q9POV3	
Isotype: IgG1	Full Name: SH3-domain binding protein 4	
Immunogen Catalog Number: AG31238	Calculated MW: 963 aa, 108 kDa	
	Observed MW: 110 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various cell lysates were subjected to SDS PAGE followed by western blot with 68045-1-Ig (SH3BP4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with Lamin B1 Monoclonal antibody (66095-1-Ig) as loading control. This data was developed using the same antibody clone with 68045-1-PBS in a different storage buffer formulation.