

For Research Use Only

TPPP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68048-5-PBS



Basic Information

Catalog Number:

68048-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19046

GenBank Accession Number:

BC131506

GeneID (NCBI):

11076

UNIPROT ID:

O94811

Full Name:

tubulin polymerization promoting protein

Calculated MW:

219 aa, 24 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G5A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50461-2, TPPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68048-4-PBS. Detection antibody: 68048-5-PBS. Standard:Ag19046. Range: 0.098-100 ng/mL.