

For Research Use Only

AIFM2/ FSP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68049-2-PBS



Basic Information

Catalog Number:

68049-2-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG30516

GenBank Accession Number:

BC023601

GeneID (NCBI):

84883

UNIPROT ID:

Q9BRQ8

Full Name:

apoptosis-inducing factor,
mitochondrion-associated, 2

Calculated MW:

41 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

4F12D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51527-1, AIFM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68049-1-PBS. Detection antibody: 68049-2-PBS.



Cytometric bead array standard curve of MP51527-1, A1FM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68049-1-PBS. Detection antibody: 68049-2-PBS. Standard: Ag30516. Range: 1.563-200 ng/mL.