

For Research Use Only

HSPB11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68059-2-PBS



Basic Information

Catalog Number:

68059-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8804

GenBank Accession Number:

BC005245

GeneID (NCBI):

51668

UNIPROT ID:

Q9Y547

Full Name:

heat shock protein family B (small),
member 11

Calculated MW:

144 aa, 16 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D5A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51201-1, HSPB11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68059-2-PBS. Detection antibody: 68059-3-PBS. Standard:Ag8804. Range: 0.098-25 ng/mL