

For Research Use Only

PYCRL Monoclonal antibody, PBS Only



Catalog Number: 68087-1-PBS

Basic Information

Catalog Number:

68087-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31938

GenBank Accession Number:

BC007993

GeneID (NCBI):

65263

UNIPROT ID:

Q53H96

Full Name:

pyrroline-5-carboxylate reductase-like

Calculated MW:

274 aa, 29 kDa

Observed MW:

26-29 kDa

Purification Method:

Protein G purification

CloneNo.:

1G7A3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Rabbit

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68087-1-Ig (PYCRL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68087-1-PBS in a different storage buffer formulation.



rabbit kidney tissue were subjected to SDS PAGE followed by western blot with 68087-1-Ig (PYCRL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68087-1-PBS in a different storage buffer formulation.