

For Research Use Only

ANXA11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68089-2-PBS



Basic Information

Catalog Number: 68089-2-PBS	GenBank Accession Number: BC007564	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 311	CloneNo.: 1C5B4
Isotype: IgG1	UNIPROT ID: P50995	
Immunogen Catalog Number: AG31839	Full Name: annexin A11	
	Calculated MW: 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP1526-1, ANXA11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68089-2-PBS. Detection antibody: 68089-3-PBS.

Cytometric bead array standard curve of MP1526-1, ANXA11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68089-2-PBS. Detection antibody: 68089-3-PBS. Standard: Ag31839. Range: 3.125-400 ng/mL.