

# FBXW11 Monoclonal antibody, PBS Only

Catalog Number: 68090-1-PBS

## Basic Information

<b>Catalog Number:</b> 68090-1-PBS	<b>GenBank Accession Number:</b> BC026213	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 23291	<b>CloneNo.:</b> 1F5A10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9UKB1	
<b>Isotype:</b> IgG1	<b>Full Name:</b> F-box and WD repeat domain containing 11	
<b>Immunogen Catalog Number:</b> AG3761	<b>Calculated MW:</b> 542 aa, 61 kDa	
	<b>Observed MW:</b> 58-62 kDa	

## Applications

**Tested Applications:**  
WB, Indirect ELISA, IHC

**Species Specificity:**  
Human, Mouse, Rat

## Background Information

FBXW11 (also known as HOS or  $\beta$ -TrCP2) is a member of F-box family proteins and plays critical role in regulating the ubiquitination of phosphorylated substrates. Abnormal expression of several FBXW11 is involved in the modulation of various biological events, such as cell cycle, differentiation, migration, inflammation, and apoptosis, through targeting multiple different substrates. For instance, FBXW11 could bind to the phosphorylated I $\kappa$ B and  $\beta$ -catenin, promoting their degradation via the ubiquitin-proteasome system. In addition, FBXW11 expression is markedly increased in mouse skin tumors and promotes tumor growth by activating the NF- $\kappa$ B signaling (PMID: 33640602). FBXW11 has 3 isoforms with the molecular mass of 58-62 kDa.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68090-1-Ig (FBXW11 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68090-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 68090-1-Ig (FBXW11 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 68090-1-PBS in a different storage buffer formulation.