

For Research Use Only

NLRP3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68102-3-PBS



Basic Information

Catalog Number:

68102-3-PBS

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26289

GenBank Accession Number:

NM_001079821

GeneID (NCBI):

114548

UNIPROT ID:

Q96P20

Full Name:

NLR family, pyrin domain containing

3

Calculated MW:

118 kDa

Purification Method:

Protein G purification

CloneNo.:

2F1B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80005-1, NLRP3 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 80422-2-PBS. Detection antibody: 68102-3-PBS. Standard:Ag26289. Range: 1.563-25 ng/mL