

For Research Use Only

# BNIP3L Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68118-4-PBS

## Basic Information

Catalog Number:

68118-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31547

GenBank Accession Number:

BC001559

GeneID (NCBI):

665

UNIPROT ID:

O60238

Full Name:

BCL2/adenovirus E1B 19kDa interacting protein 3-like

Calculated MW:

24 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D11A6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50259-1, BNIP3L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68118-2-PBS. Detection antibody: 68118-4-PBS. Standard:Ag31547. Range: 0.195-100 ng/mL.