

For Research Use Only

MAN1B1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68121-2-PBS



Basic Information

Catalog Number: 68121-2-PBS	GenBank Accession Number: BC002953	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 11253	CloneNo.: 2D6E6
Source: Mouse	UNIPROT ID: Q9UKM7	
Isotype: IgG1	Full Name: mannosidase, alpha, class 1B, member 1	
Immunogen Catalog Number: AG28302	Calculated MW: 80 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50935-1, MAN1B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68121-2-PBS. Detection antibody: 68121-3-PBS. Standard:Ag28302. Range: 0.195-100 ng/mL.