

For Research Use Only

PIN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68127-2-PBS



Basic Information

Catalog Number: 68127-2-PBS	GenBank Accession Number: BC002899	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 5300	CloneNo.: 1B7C5
Isotype: IgG1	UNIPROT ID: Q13526	
Immunogen Catalog Number: AG0767	Full Name: peptidylprolyl cis/trans isomerase, NIMA-interacting 1	
	Calculated MW: 18 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

