

For Research Use Only

# P30 of ASFV Monoclonal antibody

Catalog Number: 68134-1-Ig



## Basic Information

Catalog Number:

68134-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31446

GenBank Accession Number:

M96354

GeneID (NCBI):

22220322

UNIPROT ID:

P34204

Full Name:

Phosphoprotein p30

Calculated MW:

30 kDa

Observed MW:

35 kDa

Purification Method:

Protein G purification

CloneNo.:

3B9E3

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : Ag31446,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Recombinant protein of P30 of ASFV