

For Research Use Only

# DOHH Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68151-3-PBS



## Basic Information

Catalog Number:

68151-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13041

GenBank Accession Number:

BC002817

GeneID (NCBI):

83475

UNIPROT ID:

Q9BU89

Full Name:

deoxyhypusine  
hydroxylase/monooxygenase

Calculated MW:

33 kDa

Purification Method:

Protein G purification

CloneNo.:

2H3G4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


Selected Validation Data




Cytometric bead array sample test of MP51873-1, DOHH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68151-2-PBS. Detection antibody: 68151-3-PBS.



Cytometric bead array sample test of MP51873-2, DOHH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68151-1-PBS. Detection antibody: 68151-3-PBS.



Cytometric bead array standard curve of MP51873-1, DOHH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68151-2-PBS. Detection antibody: 68151-3-PBS. Standard:Ag13041. Range: 0.391-50 ng/mL.



Cytometric bead array standard curve of MP51873-2, DOHH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68151-1-PBS. Detection antibody: 68151-3-PBS. Standard:Ag13041. Range: 0.391-50 ng/mL.