

For Research Use Only

CACYBP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68161-4-PBS



Basic Information

Catalog Number:

68161-4-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG2340

GenBank Accession Number:

BC022352

GeneID (NCBI):

27101

UNIPROT ID:

Q9HB71

Full Name:

calcyclin binding protein

Calculated MW:

228 aa, 26 kDa

Purification Method:

Protein A purification

CloneNo.:

3D5F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51795-2, CACYBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68161-1-PBS. Detection antibody: 68161-4-PBS.



Cytometric bead array standard curve of MP51795-2, CACYBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68161-1-PBS. Detection antibody: 68161-4-PBS. Standard:Ag2340. Range: 0.781-100 ng/mL.