

For Research Use Only

# DFFA/DFF45 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68177-2-PBS



## Basic Information

Catalog Number:

68177-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG32657

GenBank Accession Number:

BC000037

GeneID (NCBI):

1676

UNIPROT ID:

O00273

Full Name:

DNA fragmentation factor, 45kDa,  
alpha polypeptide

Calculated MW:

37 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F7B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50354-1, DFFA,DFF45 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68177-1-PBS. Detection antibody: 68177-2-PBS. Standard:Ag32657. Range: 3.125-200 ng/mL