

For Research Use Only

PSMD7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68188-2-PBS



Basic Information

Catalog Number:

68188-2-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9305

GenBank Accession Number:

BC012606

GeneID (NCBI):

5713

UNIPROT ID:

P51665

Full Name:

proteasome (prosome, macropain)

26S subunit, non-ATPase, 7

Calculated MW:

324 aa, 37 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C9C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51841-1, PSMD7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68188-1-PBS. Detection antibody: 68188-2-PBS.



Cytometric bead array standard curve of MP51841-1, PSMD7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 6B188-1-PBS. Detection antibody: 6B188-2-PBS. Standard: Ag9305. Range: 0.391-50 ng/mL.