

For Research Use Only

CYP2S1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68198-6-PBS



Basic Information

Catalog Number:

68198-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4605

GenBank Accession Number:

BC033691

GeneID (NCBI):

29785

UNIPROT ID:

Q96SQ9

Full Name:

cytochrome P450, family 2, subfamily S, polypeptide 1

Calculated MW:

504 aa, 56 kDa

Purification Method:

Protein G purification

CloneNo.:

3C3C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50889-4, CYP2S1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68198-1-PBS. Detection antibody: 68198-6-PBS. Standard:Ag4605. Range: 0.391-100 ng/mL