

For Research Use Only

IDH3B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68199-2-PBS



Basic Information

Catalog Number:

68199-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8693

GenBank Accession Number:

BC001960

GeneID (NCBI):

3420

UNIPROT ID:

O43837

Full Name:

isocitrate dehydrogenase 3 (NAD+)

beta

Calculated MW:

42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C11F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50690-1, IDH3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68199-2-PBS. Detection antibody: 68199-3-PBS. Standard:Ag8693. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP50690-2, IDH3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68199-2-PBS. Detection antibody: 68199-4-PBS. Standard:Ag8693. Range: 3.125-200 ng/mL