

For Research Use Only

# IDH3B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68199-3-PBS



## Basic Information

Catalog Number:

68199-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8693

GenBank Accession Number:

BC001960

GeneID (NCBI):

3420

UNIPROT ID:

O43837

Full Name:

isocitrate dehydrogenase 3 (NAD+)

beta

Calculated MW:

42 kDa

Purification Method:

Protein G purification

CloneNo.:

1A7C8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50690-1, IDH3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68199-2-PBS. Detection antibody: 68199-3-PBS. Standard:Ag8693. Range: 1.563-200 ng/mL.