

For Research Use Only

CHCHD3/MIC19 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68259-2-PBS



Basic Information

Catalog Number:

68259-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22628

GenBank Accession Number:

BC014839

GeneID (NCBI):

54927

UNIPROT ID:

Q9NX63

Full Name:

coiled-coil-helix-coiled-coil-helix
domain containing 3

Calculated MW:

227 aa, 26 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E5B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50249-1, CHCHD3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68259-2-PBS. Detection antibody: 68259-3-PBS. Standard:Ag22628. Range: 1.563-200 ng/mL.