

For Research Use Only

FUS/TLS Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68262-2-PBS



Basic Information

Catalog Number: 68262-2-PBS	GenBank Accession Number: BC026062	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 2521	CloneNo.: 2E7D5
Isotype: IgG1	UNIPROT ID: P35637	
Immunogen Catalog Number: AG2150	Full Name: fusion (involved in t(12;16) in malignant liposarcoma)	
	Calculated MW: 75 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

