

For Research Use Only

ZHX2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68268-2-PBS



Basic Information

Catalog Number:

68268-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14420

GenBank Accession Number:

BC042145

GeneID (NCBI):

22882

UNIPROT ID:

Q9Y6X8

Full Name:

zinc fingers and homeoboxes 2

Calculated MW:

837 aa, 92 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H9E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51830-1, ZHX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68268-1-PBS. Detection antibody: 68268-2-PBS.



Cytometric bead array standard curve of MP51830-1, ZHX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 6B268-1-PBS. Detection antibody: 6B268-2-PBS. Standard: Ag14420. Range: 0.781-100 ng/mL.