

For Research Use Only

DHFR Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68274-4-PBS



Basic Information

Catalog Number:

68274-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7318

GenBank Accession Number:

BC000192

GeneID (NCBI):

1719

UNIPROT ID:

P00374

Full Name:

dihydrofolate reductase

Calculated MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

3E8E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50686-2, DHFR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68274-4-PBS. Detection antibody: 68274-5-PBS. Standard:Ag7318. Range: 0.098-50 ng/mL