

For Research Use Only

# RABEP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68285-3-PBS



## Basic Information

Catalog Number:

68285-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6023

GenBank Accession Number:

BC041700

GeneID (NCBI):

9135

UNIPROT ID:

Q15276

Full Name:

rabaptin, RAB GTPase binding effector protein 1

Calculated MW:

99 kDa

Purification Method:

Protein G purification

CloneNo.:

2D4E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

Selected Validation Data



Cytometric bead array sample test of MP51796-1, RABEP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68285-2-PBS. Detection antibody: 68285-3-PBS.



Cytometric bead array standard curve of MP51796-1, RABEP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68285-2-PBS. Detection antibody: 68285-3-PBS. Standard:Ag6023. Range: 0.781-100 ng/mL.