

# Raptor Monoclonal antibody, PBS Only

Catalog Number: 68286-1-PBS

## Basic Information

Catalog Number: 68286-1-PBS	GenBank Accession Number: BC136654	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57521	CloneNo.: 1B8A12
Source: Mouse	UNIPROT ID: Q8N122	
Isotype: IgG2b	Full Name: raptor	
Immunogen Catalog Number: AG30131	Calculated MW: 1335 aa, 149 kDa	
	Observed MW: 130-150 kDa	

## Applications

Tested Applications:  
WB, IP, Indirect ELISA

Species Specificity:  
Human

## Background Information

RPTOR, also named as KIAA1303 and RAPTOR Belongs to the WD repeat RAPTOR family. It is involved in the control of the mammalian target of rapamycin complex 1 (mTORC1) activity which regulates cell growth and survival, and autophagy in response to nutrient and hormonal signals; functions as a scaffold for recruiting mTORC1 substrates. mTORC1 is activated in response to growth factors or amino-acids. Amino-acid-signaling to mTORC1 is mediated by Rag GTPases, which cause amino-acid-induced relocalization of mTOR within the endomembrane system. Activated mTORC1 up-regulates protein synthesis by phosphorylating key regulators of mRNA translation and ribosome synthesis. mTORC1 phosphorylates EIF4EBP1 and releases it from inhibiting the elongation initiation factor 4E (eIF4E). mTORC1 phosphorylates and activates S6K1 at 'Thr-389', which then promotes protein synthesis by phosphorylating PDCD4 and targeting it for degradation.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only


For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68286-1-Ig (Raptor antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68286-1-PBS in a different storage buffer formulation.



IP result of anti-Raptor (IP:68286-1-Ig, 4ug; Detection:68286-1-Ig 1:1000) with HeLa cells lysate 1800 ug. This data was developed using the same antibody clone with 68286-1-PBS in a different storage buffer formulation.