

For Research Use Only

# KRT6A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68294-3-PBS



## Basic Information

Catalog Number: 68294-3-PBS	GenBank Accession Number: BC008807	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 3853	CloneNo.: 2G10B7
Source: Mouse	UNIPROT ID: P02538	
Isotype: IgG1	Full Name: keratin 6A	
Immunogen Catalog Number: AG0882	Calculated MW: 60 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50382-2, KRT6A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68294-3-PBS. Detection antibody: 68294-4-PBS. Standard:Ag0882. Range: 0.098-100 ng/mL.