

For Research Use Only

KRT6A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68294-4-PBS



Basic Information

Catalog Number:

68294-4-PBS

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0882

GenBank Accession Number:

BC008807

GeneID (NCBI):

3853

UNIPROT ID:

P02538

Full Name:

keratin 6A

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

1G7F5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50382-2, KRT6A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68294-3-PBS. Detection antibody: 68294-4-PBS. Standard:Ag0882. Range: 0.098-100 ng/mL.