

For Research Use Only

cTnT Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68300-5-PBS



Basic Information

Catalog Number:

68300-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7716

GenBank Accession Number:

BC002653

GeneID (NCBI):

7139

UNIPROT ID:

P45379

Full Name:

troponin T type 2 (cardiac)

Calculated MW:

36 kDa

Purification Method:

Protein G purification

CloneNo.:

3B10D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50338-3, cTnT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68300-3-PBS. Detection antibody: 68300-5-PBS. Standard:Ag7716. Range: 0.195-100 ng/mL.