

For Research Use Only

DLAT Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68303-2-PBS



Basic Information

Catalog Number: 68303-2-PBS	GenBank Accession Number: BC039084	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1737	CloneNo.: 1G6H6
Source: Mouse	UNIPROT ID: P10515	
Isotype: IgG1	Full Name: dihydrolipoamide S-acetyltransferase	
Immunogen Catalog Number: AG4226	Calculated MW: 647 aa, 70 kDa	

Applications

Tested Applications:
IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50549-1, DLAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68303-2-PBS. Detection antibody: 68303-3-PBS. Standard:Ag4226. Range: 0.098-50 ng/mL