

For Research Use Only

GTSE1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68313-3-PBS

Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 68313-3-PBS | GenBank Accession Number: BC006325 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 51512 | CloneNo.: 4B10A11 |
| Source: Mouse | UNIPROT ID: Q9NYZ3 | |
| Isotype: IgG1 | Full Name: G-2 and S-phase expressed 1 | |
| Immunogen Catalog Number: AG32817 | Calculated MW: 77 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50191-1, GTSE1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68313-2-PBS. Detection antibody: 68313-3-PBS. Standard:Ag32817. Range: 0.098-50 ng/mL