

For Research Use Only

DDT Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68314-2-PBS



Basic Information

Catalog Number: 68314-2-PBS	GenBank Accession Number: BC015508	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1652	CloneNo.: 1B7H1
Source: Mouse	UNIPROT ID: P30046	
Isotype: IgG1	Full Name: D-dopachrome tautomerase	
Immunogen Catalog Number: AG3065	Calculated MW: 118 aa, 13 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50531-1, DDT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68314-2-PBS. Detection antibody: 68314-3-PBS. Standard: Ag3065. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50531-3, DDT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68314-2-PBS. Detection antibody: 68314-4-PBS. Standard: Ag3065. Range: 0.781-100 ng/mL.