

For Research Use Only

# DDT Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68314-4-PBS



## Basic Information

Catalog Number:

68314-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3065

GenBank Accession Number:

BC015508

GeneID (NCBI):

1652

UNIPROT ID:

P30046

Full Name:

D-dopachrome tautomerase

Calculated MW:

118 aa, 13 kDa

Purification Method:

Protein G purification

CloneNo.:

3A6E8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50531-2, DDT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68314-5-PBS. Detection antibody: 68314-4-PBS. Standard:Ag3065. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50531-3, DDT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68314-2-PBS. Detection antibody: 68314-4-PBS. Standard:Ag3065. Range: 0.781-100 ng/mL