

For Research Use Only

Collagen Type III Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68320-2-PBS



Basic Information

Catalog Number:

68320-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4229

GenBank Accession Number:

BC028178

GeneID (NCBI):

1281

UNIPROT ID:

P02461

Full Name:

collagen, type III, alpha 1

Calculated MW:

1466 aa, 139 kDa

Purification Method:

Protein G purification

CloneNo.:

2A7H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50669-1, COL3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68320-2-PBS. Detection antibody: 68320-3-PBS. Standard:Ag4229. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50669-2, COL3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68320-2-PBS. Detection antibody: 68320-4-PBS. Standard:Ag4229. Range: 3.125-100 ng/mL