

For Research Use Only

EPS8L2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68330-2-PBS

Basic Information

Catalog Number:

68330-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16808

GenBank Accession Number:

BC002474

GeneID (NCBI):

64787

UNIPROT ID:

Q9H6S3

Full Name:

EPS8-like 2

Calculated MW:

715 aa, 81 kDa

Purification Method:

Protein A purification

CloneNo.:

3B10G8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50248-1, EPSL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68330-2-PBS. Detection antibody: 68330-3-PBS. Standard:Ag16808. Range: 0.781-200 ng/mL.



Cytometric bead array standard curve of MP50248-2, EPSL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68330-2-PBS. Detection antibody: 68330-4-PBS. Standard:Ag16808. Range: 0.391-200 ng/mL.