

For Research Use Only

EPS8L2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68330-4-PBS

Basic Information

Catalog Number: 68330-4-PBS	GenBank Accession Number: BC002474	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 64787	CloneNo.: 1E11G4
Source: Mouse	UNIPROT ID: Q9H6S3	
Isotype: IgG1	Full Name: EPS8-like 2	
Immunogen Catalog Number: AG16808	Calculated MW: 715 aa, 81 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50248-2, EPS8L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68330-2-PBS. Detection antibody: 68330-4-PBS. Standard:Ag16808. Range: 0.391-200 ng/mL.