

For Research Use Only

CC2D1A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68334-2-PBS



Basic Information

Catalog Number:	68334-2-PBS	GenBank Accession Number:	BC064981	Purification Method:	Protein G purification
Source:	Mouse	GeneID (NCBI):	54862	CloneNo.:	1H9F2
Isotype:	IgG1	UNIPROT ID:	Q6P1N0		
Immunogen Catalog Number:	AG10363	Full Name:	coiled-coil and C2 domain containing 1A		
		Calculated MW:	951 aa, 104 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51793-1, CC2D1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68334-1-PBS. Detection antibody: 68334-2-PBS.



Cytometric bead array standard curve of MP51793-1, CC2D1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68334-1-PBS. Detection antibody: 68334-2-PBS. Standard: Ag10363. Range: 0.781-100 ng/mL.