

For Research Use Only

RABGAP1L Monoclonal antibody

Catalog Number: 68356-1-Ig



Basic Information

Catalog Number: 68356-1-Ig	GenBank Accession Number: BC041888	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 9910	CloneNo.: 4A3A9
Source: Mouse	UNIPROT ID: B7ZAP0	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG2b	Full Name: RAB GTPase activating protein 1-like	
Immunogen Catalog Number: AG5305	Calculated MW: 815 aa, 93 kDa	
	Observed MW: 100 kDa, 140 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : U2OS cells, LNCaP cells, HeLa cells, HEK-293 cells, Jurkat cells, K-562 cells
Species Specificity: Human	

Background Information

RABGAP1L, also known as TBC1D18, is a member of the Tre2/Bub2/Cdc16 (TBC)-domain-containing protein family. Overexpression of RABGAP1L potentiated the antiviral action of IFN against a panel of viruses, including IAV, and its function required a catalytically active TBC domain.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68356-1-Ig (RABGAP1L antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.