

For Research Use Only

NDUFB7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68362-4-PBS



Basic Information

Catalog Number: 68362-4-PBS	GenBank Accession Number: BC002595	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4713	CloneNo.: 1C7E8
Source: Mouse	UNIPROT ID: P17568	
Isotype: IgG1	Full Name: NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 7, 18kDa	
Immunogen Catalog Number: AG6720	Calculated MW: 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50699-2, NDUFB7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68362-2-PBS. Detection antibody: 68362-4-PBS. Standard:Ag6720. Range: 0.195-100 ng/mL.