

For Research Use Only

CREB3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68379-2-PBS



Basic Information

Catalog Number:

68379-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1823

GenBank Accession Number:

BC009402

GeneID (NCBI):

10488

UNIPROT ID:

O43889

Full Name:

cAMP responsive element binding protein 3

Calculated MW:

44 kDa

Purification Method:

Protein G purification

CloneNo.:

2F4H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50828-1, CREB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68379-2-PBS. Detection antibody: 68379-3-PBS. Standard:Ag1823. Range: 0.098-25 ng/mL